

General Lab Product



Just Spin

Centrifuge 5420 – designed with your molecular biology routine applications in mind

This 24-place centrifuge was designed for everyday lab work. The spectrum of methods includes DNA/ RNA purification, DNA extraction from agarose gels and PCR reaction set-up. Centrifuge 5420 has a particularly small footprint to save space on the lab bench and an easy to clean keypad control panel.

You can spin up to 24 tubes from 1.5 to 2.0 mL. Optional adapters extend the volume range to 0.4 and 0.5 mL tubes as well as 0.6 mL BD Microtainer®. A rotor for 32 PCR tubes (0.2 mL) or four PCR strips is available as well.

Features and benefits

- > Automatic lid opening at the end of the run to prevent sample warming and to allow easy access to samples
- > One push short-spin key for quick centrifugation – no need to hold key all the while
- > Standard rotor provides aerosol-tight containment for the centrifugation of hazardous samples
- > Outstanding software features with menu control including 3 program buttons for your most commonly used setting
- > Particularly small footprint to save precious space on the lab bench

www.eppendorf.com/Centrifuge5420



Specifications

- > Max. RCF: 21,300 × g
- > Max. speed: 15,060 rpm
- > Max. capacity: 24 × 1.5/2.0 mL
- > Rotors option: 2 (low weight aluminum)
 - Rotor FA-24x2: 24 × 1.5/2.0 mL tubes (aerosol-tight)
 - Rotor F-32x0.2-PCR: 32 × PCR tubes (0.2 mL) or 4 × PCR strips
- > SOFT brake function: 10 ramps
- > Noise level: < 56 dB(A) with Rotor FA-24x2
- > Dimensions (W × D × H): 24 × 35 × 22 cm
- > Weight w/o accessories: 13 kg
- > Max. power consumption: 290 W

Centrifuge 5420 offers two rotor options: A 24-place standard rotor for up to 24 × 1.5 or 2.0 mL tubes and a PCR-tube rotor for 4 × PCR strip or 32 × PCR tubes (0.2 mL).

Ordering Information

Centrifuges	Order no.
Centrifuge 5420, non-refrigerated, without rotor, keypad, 230 V/50–60 Hz	5420 000 318
Centrifuge 5420, non-refrigerated, with aerosol-tight Rotor FA-24x2, keypad, 230 V/50–60 Hz	5420 000 113
Rotors & Adapters	
Rotor FA-24x2 for Centrifuge 5420, aerosol-tight, for 24 × 1.5/2.0 mL tubes, incl. aerosol-tight rotor lid	5495 400 001
Rotor F-32x0.2-PCR for Centrifuge 5420, for 32 × 0.2 mL PCR tubes or 4 × 8 PCR strips, incl. rotor lid	5495 404 007
Adapter, for 1 tube 0.5 mL or microtainer, max. Ø 8 mm, for all 1.5/2.0 mL rotors, set of 6 pcs.	5425 716 001
Adapter, for 1 PCR tube 0.2 mL, max. Ø 6 mm, for all 1.5/2.0 mL rotors, set of 6 pcs.	5425 715 005
Adapter, for 1 tube 0.4 mL, max. Ø 6 mm, for all 1.5/2.0 mL rotors, set of 6 pcs.	5425 717 008

S.I.A.L. Srl
Via Giovanni Devoti 14
00167 - Roma, Italia

Tel: +39 06/6625280
Fax: +39 06/6628503
E-mail: info@sialgroup.com

Orari: Lun - Ven 8:30 - 18:00
Web: sialgroup.com

CHIAMACI AL NUMERO
06 6625280