

Cole-Parmer®

 **Sial**

Cole-Parmer® Colony Counter CC-200

- Pressure sensitive counting
- Average count facility
- Bright white energy saving LED lighting
- With BioCote antimicrobial protection
- Audible count confirmation
- Choice of light or dark background
- Connectivity to printer or computer





Cole-Parmer® Colony Counter CC-200

Touch pressure with felt tip marker on petri dish registers cumulative count on the digital display with confirmation by audible tone (can be turned on or off). The pressure required to register a count can be adjusted to suit each user.

Averaging facility calculates average count over multiple plates. Counting results as well as useful statistics including SD can be sent directly to the accessory printer or to a computer via a USB cable supplied.

Sub-stage illumination by low energy bright LED's allows glare-free optimum viewing. A switchable black background is provided to enhance viewing of translucent and difficult to see colonies. Supplied with two Wolffhuegel graticules and dish centering adaptors to facilitate use with 50mm to 90mm dishes. A choice of magnifiers and a printer are available as optional accessories.

Technical Specification

Specification	CC-200
Lighting	White LED array
Digital display	3 digit LED
Count	0 to 999
Dimensions (w x d x h)	310 x 300 x 140mm
Net weight	1.5kg
Electrical supply	90 to 240V, 50/ 60Hz, 70W

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Colony counter, advanced	14212-04	CC-200	CC-200	SC6PLUS
Colony counter, advanced with IQ/OQ documentation	14212-06	CC-200	CC-200	SC6PLUS/IQOQ
1.7x magnifier	14212-60	CC-200	N/A	SC6/ 1
3x magnifier	14212-61	CC-200	N/A	SC6/ 1/ 3
Wolffhuegel graticule & segmentation discs (pack of 10)	99952-84	CC-200	N/A	SC6/ 2
Spare dish centering adaptors (pack of 2)	99952-85	CC-200	N/A	SC6/ 3
Clear protective discs (pack of 5)	99952-86	CC-200	N/A	SC6/ 4
Spare receiving dish	99952-87	CC-200	N/A	SC6/ 5
Accessory Printer	03012-99	CC-200	N/A	SMP30/ 1

coleparmer.com

Sial

S.I.A.L. Srl
Via Giovanni Devoti 14
00167 - Roma, Italia

Tel: +39 06/6625280
Fax: +39 06/6628503
E-mail: info@sialgroup.com

Orari: Lun - Ven 8:30 - 18:00
Web: sialgroup.com

CHIAMACI AL NUMERO
06 6625280

Cole-Parmer®

essentials

an antylia scientific company

Find out more at antylia.com